



## Montrose Specifications

<b>Montrose Major</b>		
<b>Seat Angle</b>	<b>Flat Seat Height (floor to top of seat)</b>	<b>Ramped Seat Height (floor to top of seat)</b>
67°	535mm	565mm
74°	495mm	525mm
81°	455mm	485mm
88°	415mm	445mm

<b>Montrose Minor</b>		
<b>Seat Angle</b>	<b>Flat Seat Height (floor to top of seat)</b>	<b>Ramped Seat Height (floor to top of seat)</b>
67°	485mm	515mm
74°	445mm	475mm
81°	405mm	435mm
88°	365mm	395mm

<b>Dimensions (mm)</b>	<b>Montrose Major</b>	<b>Montrose Minor</b>
Seat Width	470mm	420mm
Seat Length	380-500mm	330-450mm
Back Recline	5°-32°	5°-32°
Back Height (from seat)	780mm	730mm
Arm Height (from seat)	240mm	240mm
Overall Width	680mm	630mm
Overall Length	950mm	900mm
Overall Height	1190mm	1090mm
Overall Arm Height (from floor)	680mm	630mm
<b>Legrest</b>		
Height	455mm	405mm
Width	580mm	530mm
Length	480mm	480mm