

Kirton H20 Shower, Toilet Commode Chair Specification Guide

	Neptune Shower Stool	M2 Shower, Toilet & Commode Chair	M2 Self-Propelled Shower, Toilet & Commode Chair	M2 Gas-Tip Shower, Toilet & Commode Chair	M2 EI-Tip Shower, Toilet & Commode Chair	M2 Mini EI-Tip Shower, Toilet & Commode Chair	200kg Shower Stool	M2 200kg Shower, Toilet & Commode Chair	M2 Motorised 200kg Shower, Toilet & Commode Chair
Product Code	HM001	HM003	HM004	HM007	HM006	HM009	HM002	HM008	HM005
Weight Limit	23 stones 150kg	20 stones 125kg	20 stones 125kg	22 stones 140kg	22 stones 140kg	12 stones 80kg	30 stones 200kg	30 stones 200kg	30 stones 200kg
Weighth of chair	Stool: 10lbs/4.6kg Backrest: 1.3 kg Complete: 5.9 kg.	2 stones 7lbs 16kg	2 stones 10lbs 17.6kg	3 stones 4lbs 21kg	4 stones 25.8kg	3 stones 10lbs 23.6kg	1 stone 6lbs 9.4kg	2 stones 10lbs 17.5kg	4 stones 10lbs 30kg
Specifications (mm)									
Seat width	410 500*	465	465	465	465	460	700	560	560
Seat depth	365	440	440	440	440	375	360	440	440
Seat height(s)	430-640 (8 adjustment positions)	445-565 (7 adjustment positions)	470-490-560 (3 adjustment positions)	555-615 (4 adjustment positions)	555-615 (4 adjustment positions)	460-540 (5 adjustment positions)	470-500-530 (3 adjustment positions)	565	565
Width between armrests	410	460	460	460	460	340	700	560	560
Height (without headrest)	560-770 (to armrest)	840-970	880-900-970	1000-1060	1000-1060	1260 (when handles in highest position)	610-670	970	970
Overall width	545-585 (at floor)	570 armrest to armrest .	630	570	590	595	830	680	715
Overall length (in upright position)	440-520 (at floor)	660	775	770	745	655	480	685	710
Tilt-in-Space	n/a	n/a	n/a	37 back; 7 forwards	37 back; 7 forwards	37 back; 7 forwards	n/a	n/a	n/a
Toilet height clearance	n/a	355-485	380-400-465	465-525	465-525	390-470	n/a	465	460
* available on request									
Dimensions in mm									